Clinical and Analytical Performance of a New Molecular C. difficile Direct Assay

Poster Number : 312 ASM Clinical Virology Symposium Savannah, GA I May 7-10, 2017

Revised Abstract

Introduction: Clostridium difficile, a gram-positive, anaerobic, sporeforming bacillus, is associated with significant morbidity and mortality. C. difficile has been linked to 50% to 75% of antibiotic-associated colitis cases, up to one-third of all antibiotic-associated diarrhea cases and over 90% of all cases of antibiotic-associated pseudomembranous colitis.

The DiaSorin Molecular Simplexa[®] C. difficile Direct Assay is FDAcleared for sample-to-answer detection of the *Clostridium difficile* toxin B gene (tcdB) from liquid or unformed stool samples from individuals suspected of difficile infection (CDI). This test utilizes realtime PCR amplification without any nucleic acid extraction.

The analytical performance characteristics, including analytical sensitivity/limit of detection (LoD), reproducibility, analytical reactivity, analytical specificity/cross-reactivity, microbial inhibition and interfering substances, were determined. Clinical agreement studies compared Simplexa performance to direct and combination direct/enriched culture and toxin assay. Positive and negative agreement were also determined for Simplexa versus three FDA cleared NAAT assays.

Methods: The analytical sensitivity/limit of detection was tested for two toxigenic C. difficile strains, ATCC 43255 and NAP1A. A reproducibility panel with both strains was tested at 3 different laboratories. For the method comparison study, the performance of Simplexa C. difficile Direct was compared to direct culture and direct/enriched combined culture with toxin assay. Samples from the method comparison study were also tested with one of three FDA cleared molecular NAAT assays. Analytical reactivity was determined for an additional 32 toxigenic *C. difficile* strains. In order to evaluate analytical specificity, a panel of 127 bacteria, viruses and fungi, were tested for cross-reactivity. Microbial inhibition was evaluated by testing low levels of ATCC 43255 or NAP1A in a background of high levels of each of the same 127 microorganisms. The potential inhibitory effect of 33 endogenous and exogenous substances on detection of ATCC 43255 and NAP1A was evaluated.

Results: The LoD for the Simplexa[®] C. difficile Direct Assay was 0.95 CFU/mL for the ATCC 43255 strain and 0.43 CFU/mL for NAP1A. Across all 3 sites in the reproducibility study, 100% detection was observed for ATCC 43255 LP, ATCC 43255 MP and NAP1A MP. For NAP1A LP, the detection rate was 98.9%. In comparison to direct culture and toxin assay, Simplexa C. difficile Direct had a positive agreement of 95.7% and a negative agreement of 92.1%. When Simplexa C. difficile Direct was compared to direct/enriched culture and toxin assay, the sensitivity was 85.9% and the specificity was 95.1%. The performance of Simplexa C. difficile Direct compared to 3 FDA cleared NAAT assays, demonstrated 93.4%, 93.9% and 84.8% positive agreement and 96.6%, 94.0% and 99.2% negative agreement. Simplexa C. difficile Direct was able to detect all 32 toxigenic strains tested for analytical reactivity. The assay did not have any cross-reactivity to the 127 microorganisms tested; at the same time, none of the organisms tested inhibited the detection of the ATCC 43255 or NAP1A strains when tested at 2-4x LoD. The assay was able to detect ATCC 43255 and NAP1A in the presence of each potentially interfering substance tested.

Real-time PCR Amplification and Detection: Simplexa C. difficile Direct Assay (MOL2950) contains all reagents for onboard extraction and real-time PCR. Swabs were dipped into stool specimens and excess liquid was removed. The swab was placed into Sample Prep Buffer and swirled to release the stool. For each reaction on the Direct Amplification Disc (DAD), 50 μ L of sample was loaded into the sample port and 50 μ L of C. *difficile* Direct Reaction Mix was loaded into the reaction port. All testing was performed using the LIAISON MDx instrument.

Limit of Detection: The LoD for C. difficile ATCC 43255 and NAP1A (Zeptometrix) strains was determined as the lowest concentration with $\geq 95\%$ detection in negative stool-TE matrix.

Reproducibility: The panel included a positive control (PC), a negative control and 4 contrived samples: a low positive (LP) and a medium positive (MP) for the *C. difficile* ATCC 43255 and NAP1A strains. The *C. difficile* reproducibility panel was contrived in negative stool-TE matrix.

Method Comparison: A panel of clinical specimens from 5 geographically diverse sites was evaluated using the Simplexa Direct and Direct/Enriched Culture methods as well as 3 commercially available FDA-cleared molecular Nucleic Acid Amplification Test (NAAT) assays. Positive and negative agreements were determined.

Analytical Reactivity: Thirty-two toxigenic C. difficile strains were individually spiked into negative stool-TE matrix at 2-4X LoD and tested.

Cross-reactivity: The test panel consisted of 10⁶ CFU/mL for bacteria and fungi, 10^6 cells/mL for parasites or 10^5 TCID₅₀/mL for viruses and was formulated in negative stool-TE matrix.

Microbial Inhibition: The organisms in the cross-reactivity panel were spiked into a separate baseline sample of *C. difficile* ATCC 43255 or NAP1A at 2-4X LoD in negative stool-TE matrix.

Potentially Interfering Substances: Each individual substance was spiked into a baseline sample consisting of *C. difficile* ATCC 43255 or NAP1A at 2-4X LoD in negative stool-TE matrix.

C. diff

L. Geller, C. Austria, B. Yu, S. Penilla, C. Timmerman, S. McCloud, Y. Parocua, A. Bologa and M. Tabb DiaSorin Molecular LLC, Cypress, CA, USA

Methods

Results

Limit of Detection was 0.95 CFU/mL or less for **C. difficile strains tested in Table 1.**

Table 1. Limit of Detection

C. difficile Bacterial Strain	LoD Concentration (CFU/mL)
ATCC 43255	0.95
NAP1A	0.43

Results (continued)

Reproducibility: All ATCC 43255 replicates were detected at all sites; for NAP1A, 100% detection for MP and 98.9% detection (89/90 replicates) for LP was observed. PC and negative control were as expected.

≥ 92.1% (Tables 2 – 6).

Table 2. Simplexa vs. Direct Culture

	Direct Culture + Toxin Assay				
Simplexa Results	Detected Not Detected Total				
Detected	265	163	428		
Not Detected	12	1890	1902		
Total	277	2053	2330		
Positive Agreement	95.7% (265/277)	Nogotivo Agroomont	92.1% (1890/2053)		
	95% CI: 92.6% to 97.5%	Negative Agreement	95% CI: 90.8% to 93.2%		

Table 3. Simplexa vs. Combined Culture

	Combined Culture (Direct & Enriched) + Toxin Assay		
Simplexa Results	Detected	Not Detected	Total
Detected	336	96	432
Not Detected	55	1849	1904
Total	391	1945	2336
Constitute	85.9% (336/391)	Specificity	95.1% (1849/1945)
Sensitivity	95% CI: 82.1% to 89.0%	L% to 89.0%	95% CI: 94.0% to 95.9%

Table 4. Simplexa vs. NAAT 1

	NAAT1		
Simplexa Results	Detected	Not Detected	Total
Detected	114	24	138
Not Detected	8	683	691
Total	122	707	829
Positive Agreement	93.4% (114/122)	Nagativo Agroomont	96.6% (683/707)
	95% CI: 87.6% to 96.6%	Negative Agreement	95% CI: 95.0% to 97.7%

Table 5. Simplexa vs. NAAT 2

	NAAT2		
Simplexa Results	Detected	Not Detected	Total
Detected	138	38	176
Not Detected	9	591	600
Total	147	629	776
Positive Agreement	93.9% (138/147)	Nogotivo Agroomont	94.0% (591/629)
	95% CI: 88.8% to 96.7%	Negative Agreement	95% CI: 91.8% to 95.6%

Table 6. Simplexa vs. NAAT 3

	NAAT3		
Simplexa Results	Detected	Not Detected	Total
Detected	112	5	117
Not Detected	20	584	604
Total	132	589	721
Positive Agreement	84.8% (112/132)	Nogotivo Agroomont	99.2% (584/589)
	95% CI: 77.8% to 90.0%	Negative Agreement	95% CI: 98.0% to 99.6%

Method Comparison – all comparators: Positive Agreement was \geq 84.8%. Negative Agreement was

Results (continued)

Analytical Reactivity: All toxigenic *C. difficile* strains tested were detected at 512 CFU/mL or lower.

Table 7. C. difficile Strains Tested

Strain	Toxinotype	Toxin	Ribotype
ATCC 17857	0	A+B+	001
ATCC 17858	0	A+B+	054
ATCC 43594	0	A+B+	005
ATCC 43596	0	A+B+	012
ATCC 43598	VIII	A-B+	017
ATCC 43599	0	A+B+	001
ATCC 43600	0	A+B+	014
ATCC 51695	0	A+B+	001
ATCC 700792	0	A+B+	005
ATCC 9689	0	A+B+	001
BAA-1382	0	A+B+	012
BAA-1805	IIIb	A+B+	027
BAA-1814	XXII	A+B+	251
BAA-1870	IIIb	A+B+	027
BAA-1871	0	A+B+	001
BAA-1872	0	A+B+	207
BAA-1873	0	A+B+	053
BAA-1874	0	A+B+	002
BAA-1875	V	A+B+	078
CCUG 8864 (CCUG 20309)	X	A-B+	Unknown
IS81		A+B+	034
R1880		A+B+	086
R7771	VIII	A-B+	110
R8366	0	A+B+	001
R8637	IX	A+B+	019
R9367	XIII	A+B+	070
R9385	XV	A+B+	122
R10456	IV	A+B+	058
R10725	V	A+B+	078
R10842	VI	A+B+	045
R10870	XIV	A+B+	111
R12425		A+B+	103

Refer to the Simplexa *C. difficile* Direct Package Insert for information on the concentration tested for each strain.

No interference was detected with the substances tested in Table 8.

Table 8. Potentially Interfering Substances Tested

Afrin	Dramamine	Metronidazole	Palmitic Acid
Antacid and Anti-gas	Dulcolax	Milk of Magnesia	Pepto-Bismol
Anusol Plus	Gynol II, Nonoxynol 9	Mineral Oil	Peparation H
Barium Sulfate	Hydrocortisone cream	Monistat 7	Sennosides
Benzalkonium	Imodium	Mucin	SPF 30 Sunscreen
Blood	KY Jelly	Naproxen	SPF 50 Sunscreen
Colace	Mesalazine	Nystatin	Stearic Acid

Refer to the Simplexa C. difficile Direct Package Insert for information on the concentration tested for each strain.

Correspondence: Louis Geller Louis.Geller@diasorin.com

Tums Vancomy Vaseline Witch haze

Results (continued)

No cross-reactivity or microbial inhibition was observed with the organisms tested in Table 9

Table 9. Cross Reactivity Pathogens (representative of the 127 organisms tested).

Adenovirus 1	Clostridium sepiticum	Lactobacillus acidophilus
Alcaligenes faecalis subsp. faecalis	Clostridium sordellii	Listeria monocytogenes
Bacillus cereus	Clostridium sporogenes	Norovirus G2
Bacteroides fragilis	Clostridium tetani	Proteus mirabilis
Campylobacter coli	Coxsackievirus A10	Pseudomonas aeruginosa
Campylobacter jejuni	Cytomegalovirus	Rotavirus
Candida albicans	Echovirus 11	<i>Salmonella enterica</i> subsp. arizonae
Chlamydia trachomatis	Enterobacter aerogenes	<i>Salmonella enterica</i> subsp. enterica
Citrobacter freundii	Enterobacter cloacae	Serratia liquefaciens
Citrobacter koseri	<i>Enterococcus faecalis</i> vanB	Shigella sonnei
<i>Clostridium difficile</i> (non-toxigenic)	<i>Enterococcus faecium</i> vanA	Staphylococcus aureus
Clostridium haemolyticum	Enterovirus 71	Staphylococcus epidermidis
Clostridium innocuum	Escherichia coli	Streptococcus agalactiae
Clostridium perfringens	Helicobacter pylori	Vibrio cholerae

Refer to the Simplexa *C. difficile* Direct Package Insert for the complete list of organisms tested for cross reactivity and microbial inhibition.

Conclusions

- The Simplexa C. difficile Direct kit and LIAISON[®] MDx from DiaSorin Molecular provide sample-toanswer detection of *Clostridium difficile* from unextracted liquid or unformed stool in approximately one hour, with a simple workflow.
- Simplexa *C. difficile* Direct does not show crossreactivity to, nor inhibition by, other organisms found in human stool. The assay could detect all toxigenic *C. difficile* strains tested.
- This assay is FDA-cleared with a CLIA moderate complexity rating and is CE-marked.



www.diasorin.com www.focusdx.com