



Evaluation of Two Molecular Assays for the Detection of *Clostridium difficile*

Matthew England, Patrick Erdman, and Wallace Greene

Diagnostic Virology Laboratory, Department of Pathology and Laboratory Medicine
Penn State Health Milton S. Hershey Medical Center, Hershey, PA

Penn State Health
Milton S. Hershey
Medical Center
717-531-5121
whg2@psu.edu

Introduction

Clostridium difficile causes potentially serious cases of infectious diarrhea and is a common hospital acquired infection (HAI), with nearly half a million individuals being infected each year. Recently, several nucleic acid amplification techniques have been developed to improve sensitivity for detecting *C. difficile*. We compared the Meridian *illumigene*® *C. difficile* and the recently FDA-cleared Focus Simplexa™ *C. difficile* Direct assays to standard culture methods. The *illumigene*® assay uses loop-mediated isothermal DNA amplification (LAMP) targeting the PaLoc gene locus, whereas Simplexa™ uses real-time PCR targeting the toxin B gene. Both assays were found to detect more *C. difficile* than toxigenic culture; however, the Simplexa™ Direct assay was more sensitive than culture and *illumigene*®.

Methods

- In this study, 410 unformed stool specimens were received in the laboratory and tested prospectively with the *illumigene*® and Simplexa™ *C. difficile* assays according to manufacturers' instructions.
- Aliquots of all specimens were sent to a reference laboratory for toxigenic culture by standard methods.
- Results of the molecular assays were compared to each other as well as to culture.

illumigene® *C. difficile* Assay Test Procedure



Simplexa™ *C. difficile* Direct Assay Test Procedure



Results

Performance of *illumigene*® and Simplexa™ *C. difficile* Assays

	Simplexa™ Positive	Simplexa™ Negative
<i>illumigene</i> ® Positive	67 ^a	2 ^c
<i>illumigene</i> ® Negative	13 ^b	328 ^d

^a – 9 were negative by culture

^b – 11 were culture positive, 2 were culture negative

^c – Both were culture negative

^d – 3 were culture positive

Performance of Simplexa™ *C. difficile* Assay vs Culture

	Simplexa™	
Culture	Positive	Negative
Positive	69	3
Negative	11	327

Performance of *illumigene*® *C. difficile* Assay vs Culture

	<i>illumigene</i> ®	
Culture	Positive	Negative
Positive	58	14
Negative	11	327

Sensitivity and Specificity^a

	<i>illumigene</i> ®	Simplexa™
Sensitivity	81 %	96 %
Specificity	97 %	97 %

^a – with culture as the gold-standard

Conclusions

- The two molecular assays were in agreement 395/410 or 96 % of the time.
- The molecular assays detected 9 positive samples that were negative by culture.
- Culture detected 3 positive samples that were negative by both molecular assays.
- The Simplexa™ assay detected 13 positive samples negative by the *illumigene*® assay – 11 confirmed by culture, 2 were not.
- The *illumigene*® assay detected 2 positive samples that were negative by the Simplexa™ assay; however both were negative by culture.
- Overall, the Simplexa™ assay was more sensitive than the *illumigene*® assay, but both were more sensitive than toxigenic culture.
- This Simplexa™ *C. difficile* Direct Assay is a sensitive and simple method to detect the presence of *C. difficile* in stool.